

# Analysis of high-throughput sequencing data


(or: how to edit Very Large Files)



# Cluster generation




1. PREPARE GENOMIC DNA SAMPLE




Randomly fragment genomic DNA and ligate adapters to both ends of the fragments.

2. ATTACH DNA TO SURFACE




Bind single-stranded fragments randomly to the inside surface of the flow cell channels.

3. BRIDGE AMPLIFICATION




Add unlabeled nucleotides and enzyme to initiate solid-phase bridge amplification.

4. FRAGMENTS BECOME DOUBLE STRANDED




The enzyme incorporates nucleotides to build double-stranded bridges on the solid-phase substrate.

5. DENATURE THE DOUBLE-STRANDED MOLECULES



Denaturation leaves single-stranded templates anchored to the substrate.


6. COMPLETE AMPLIFICATION



Several million dense clusters of double-stranded DNA are generated in each channel of the flow cell.

# Sequencing

7. DETERMINE FIRST BASE




First chemistry cycle: to initiate the first sequencing cycle, add all four labeled reversible terminators, primers and DNA polymerase enzyme to the flow cell.

8. IMAGE FIRST BASE




After laser excitation, capture the image of emitted fluorescence from each cluster on the flow cell. Record the identity of the first base for each cluster.

9. DETERMINE SECOND BASE




Second chemistry cycle: to initiate the next sequencing cycle, add all four labeled reversible terminators and enzyme to the flow cell.

10. IMAGE SECOND CHEMISTRY CYCLE




After laser excitation, collect the image data as before. Record the identity of the second base for each cluster.

11. SEQUENCE READS OVER MULTIPLE CHEMISTRY CYCLES



Repeat cycles of sequencing to determine the sequence of bases in a given fragment a single base at time.

12. ALIGN DATA



Align data, compare to a reference, and identify sequence differences.

40 million clusters per flow cell



20 microns

illumina

# Now what?

- Download FASTQ files from the sequencing center.
- VERY LARGE – do not try to open in a text editor.

[illegible]

# FASTQ files

# 1 read = 4 lines

@HWI-ST141\_0365:2:1101:2983:2114#TTAGGC/1  
TGAGAATGTGAAGGCCAAGATCCAAGACAAGGAAGGGATTCCCCAGACC  
+HWI-ST141\_0365:2:1101:2983:2114#TTAGGC/1  
Y^\_ccZz\\\b``Rbfb`Z`^`eaRbReJyBHYBH^c]YXa^\_\_^acc  
@HWI-ST141\_0365:2:1101:2965:2155#TTAGGC/1  
GGCATCTGACAGTTGTATTTGAGATGGTATGCCACACACCTAGTTAAAGTA  
+HWI-ST141\_0365:2:1101:2965:2155#TTAGGC/1  
\_bbcreeeeegghhdefihhacgggf`eggghhgiiiiihaedhhfd  
@HWI-ST141\_0365:2:1101:2886:2184#TTAGGC/1  
CTCCACATAGCTGAACGTTGATACAGATCTACAAAAATAATGAAATGAT  
+HWI-ST141\_0365:2:1101:2886:2184#TTAGGC/1  
bbbeeeeeggggiiiiiiihihihihihihihihihihihihihihi  
@HWI-ST141\_0365:2:1101:2778:2215#TTAGGC/1  
CGAAAACCACTTCTCAGGGAGGGGAGAGTGCATACATGTTGTC  
+HWI-ST141\_0365:2:1101:2778:2215#TTAGGC/1  
bbbeeeeeggggihihiifgiiehf`gghiiiiiiiiii  
@HWI-ST141\_0365:2:1101:2956:2216#TTAGGC/1  
CATGGATGCTCTCAAAGTGTGTCTGATATGGGCTACTACATCGAGGACA  
+HWI-ST141\_0365:2:1101:2956:2216#TTAGGC/1  
\_bbeeeeeggggiiifbgfghfghagffgghfghfiifcbgaaffcf  
@HWI-ST141\_0365:2:1101:2985:2249#TTAGGC/1  
GAGAAACCTCCGACACTGGCTG  
+HWI-ST141\_0365:2:1101:2985:2249#TTAGGC/1  
Bbbeeeeegggfhhhhhhhh

.....  
.....



# FASTQ files

**@HWI-ST141\_0365:2:1101:2983:2114#TTAGGC/1**

TGAGAATGTGAAGGCCAAGATCCAAGACAAGGAAGGGATTCCCCCAGACC

+HWI-ST141\_0365:2:1101:2983:2114#TTAGGC/1

Y^\_ccZZa\\b`Rbfb`Z`^`eaRbReJYbHYYbH^c]YXa^\_\_^acc

- First line is an IDENTIFIER.
- Starts with @, then instrument name, flowcell lane, tile number and flowcell x,y coordinates.
- Ends with the barcode sequence and pair # (for paired-end sequencing).

# FASTQ files

```
@HWI-ST141_0365:2:1101:2983:2114#TTAGGC/1
TGAGAATGTGAAGGCCAAGATCCAAGACAAGGAAGGGATTCCCCCAGACC
+HWI-ST141_0365:2:1101:2983:2114#TTAGGC/1
Y^_ccZZa\\b`Rbfb`Z`^`eaRbReJYbHYYbH^c]YXa^__^acc
```

- Second line is the READ SEQUENCE

# FASTQ files

```
@HWI-ST141_0365:2:1101:2983:2114#TTAGGC/1
TGAGAATGTGAAGGCCAAGATCCAAGACAAGGAAGGGATTCCCCCAGACC
+HWI-ST141_0365:2:1101:2983:2114#TTAGGC/1
Y^_ccZZa\\b`Rbfb`Z`^`eaRbReJYbHYYbH^c]YXa^__^acc
```

- Third line is the same IDENTIFIER, now starting with + (for additional information)



# FASTQ files

```
@HWI-ST141_0365:2:1101:2983:2114#TTAGGC/1
TGAGAATGTGAAGGCCAAGATCCAAGACAAGGAAGGGATTCCCCCAGACC
+HWI-ST141_0365:2:1101:2983:2114#TTAGGC/1
Y^_ccZZa\\b`Rbfb`Z`^`eaRbReJYbHYYbH^c]YXa^__^acc
```

- Fourth line is Phred scale QUALITY SCORES.
- Based on the computer ASCII numbers.
- Range from 0 to 40.
- $Q = -10 \log P$ , where  $P$  is the probability of incorrect base call.

<u>Phred Quality Score</u>	<u>Probability of incorrect base call</u>	<u>Base call accuracy</u>
10	1 in 10	90 %
20	1 in 100	99 %
30	1 in 1000	99.9 %
40	1 in 10000	99.99 %


# FASTQ files

- All keyboard characters represent an ASCII number.
- Saves space, e.g. U (1 byte) instead of 85 (2 bytes).
- For Illumina data, ASCII 33 (!) = 0, ASCII 73 (I) = 40

Dec	Hx	Oct	Char	Dec	Hx	Oct	Html	Chr	Dec	Hx	Oct	Html	Chr	Dec	Hx	Oct	Html	Chr
0	0	000	<b>NUL</b> (null)	32	20	040	&#32;	Space	64	40	100	&#64;	@	96	60	140	&#96;	`
1	1	001	<b>SOH</b> (start of heading)	33	21	041	&#33;	!	65	41	101	&#65;	A	97	61	141	&#97;	a
2	2	002	<b>STX</b> (start of text)	34	22	042	&#34;	"	66	42	102	&#66;	B	98	62	142	&#98;	b
3	3	003	<b>ETX</b> (end of text)	35	23	043	&#35;	#	67	43	103	&#67;	C	99	63	143	&#99;	c
4	4	004	<b>EOT</b> (end of transmission)	36	24	044	&#36;	\$	68	44	104	&#68;	D	100	64	144	&#100;	d
5	5	005	<b>ENQ</b> (enquiry)	37	25	045	&#37;	%	69	45	105	&#69;	E	101	65	145	&#101;	e
6	6	006	<b>ACK</b> (acknowledge)	38	26	046	&#38;	&	70	46	106	&#70;	F	102	66	146	&#102;	f
7	7	007	<b>BEL</b> (bell)	39	27	047	&#39;	'	71	47	107	&#71;	G	103	67	147	&#103;	g
8	8	010	<b>BS</b> (backspace)	40	28	050	&#40;	(	72	48	110	&#72;	H	104	68	150	&#104;	h
9	9	011	<b>TAB</b> (horizontal tab)	41	29	051	&#41;	)	73	49	111	&#73;	I	105	69	151	&#105;	i
10	A	012	<b>LF</b> (NL line feed, new line)	42	2A	052	&#42;	*	74	4A	112	&#74;	J	106	6A	152	&#106;	j
11	B	013	<b>VT</b> (vertical tab)	43	2B	053	&#43;	+	75	4B	113	&#75;	K	107	6B	153	&#107;	k
12	C	014	<b>FF</b> (NP form feed, new page)	44	2C	054	&#44;	,	76	4C	114	&#76;	L	108	6C	154	&#108;	l
13	D	015	<b>CR</b> (carriage return)	45	2D	055	&#45;	-	77	4D	115	&#77;	M	109	6D	155	&#109;	m
14	E	016	<b>SO</b> (shift out)	46	2E	056	&#46;	.	78	4E	116	&#78;	N	110	6E	156	&#110;	n
15	F	017	<b>SI</b> (shift in)	47	2F	057	&#47;	/	79	4F	117	&#79;	O	111	6F	157	&#111;	o
16	10	020	<b>DLE</b> (data link escape)	48	30	060	&#48;	0	80	50	120	&#80;	P	112	70	160	&#112;	p
17	11	021	<b>DC1</b> (device control 1)	49	31	061	&#49;	1	81	51	121	&#81;	Q	113	71	161	&#113;	q
18	12	022	<b>DC2</b> (device control 2)	50	32	062	&#50;	2	82	52	122	&#82;	R	114	72	162	&#114;	r
19	13	023	<b>DC3</b> (device control 3)	51	33	063	&#51;	3	83	53	123	&#83;	S	115	73	163	&#115;	s
20	14	024	<b>DC4</b> (device control 4)	52	34	064	&#52;	4	84	54	124	&#84;	T	116	74	164	&#116;	t
21	15	025	<b>NAK</b> (negative acknowledge)	53	35	065	&#53;	5	85	55	125	&#85;	U	117	75	165	&#117;	u
22	16	026	<b>SYN</b> (synchronous idle)	54	36	066	&#54;	6	86	56	126	&#86;	V	118	76	166	&#118;	v
23	17	027	<b>ETB</b> (end of trans. block)	55	37	067	&#55;	7	87	57	127	&#87;	W	119	77	167	&#119;	w
24	18	030	<b>CAN</b> (cancel)	56	38	070	&#56;	8	88	58	130	&#88;	X	120	78	170	&#120;	x
25	19	031	<b>EM</b> (end of medium)	57	39	071	&#57;	9	89	59	131	&#89;	Y	121	79	171	&#121;	y
26	1A	032	<b>SUB</b> (substitute)	58	3A	072	&#58;	:	90	5A	132	&#90;	Z	122	7A	172	&#122;	z
27	1B	033	<b>ESC</b> (escape)	59	3B	073	&#59;	;	91	5B	133	&#91;	[	123	7B	173	&#123;	{
28	1C	034	<b>FS</b> (file separator)	60	3C	074	&#60;	<	92	5C	134	&#92;	\	124	7C	174	&#124;	
29	1D	035	<b>GS</b> (group separator)	61	3D	075	&#61;	=	93	5D	135	&#93;	]	125	7D	175	&#125;	}
30	1E	036	<b>RS</b> (record separator)	62	3E	076	&#62;	>	94	5E	136	&#94;	^	126	7E	176	&#126;	~
31	1F	037	<b>US</b> (unit separator)	63	3F	077	&#63;	?	95	5F	137	&#95;	_	127	7F	177	&#127;	DEL

# Quality control processing of raw data

- We need to make sure that no low quality reads are in the data
- Also that no adapter sequences are left in the dataset





- We also want to study the quality score distribution (should be high) and nucleotide distribution (should be random), and the fraction of duplicate reads (should be low).



# fastQC

Quick way to visualize raw data, automatically generates html report file.



# The FASTX toolkit

FASTQ-to-FASTA converter

Convert FASTQ files to FASTA files.

FASTQ Information

Chart Quality Statistics and Nucleotide Distribution

FASTQ/A Collapser

Collapsing identical sequences in a FASTQ/A file into a single sequence (while maintaining reads counts)

FASTQ/A Trimmer

Shortening reads in a FASTQ or FASTQ files (removing barcodes or noise).

FASTQ/A Renamer

Renames the sequence identifiers in FASTQ/A file.

FASTQ/A Clipper

Removing sequencing adapters / linkers

# The FASTX toolkit

## FASTQ/A Reverse-Complement

Producing the Reverse-complement of each sequence in a FASTQ/FASTA file.

## FASTQ/A Barcode splitter

Splitting a FASTQ/FASTA files containing multiple samples

## FASTA Formatter

changes the width of sequences line in a FASTA file

## FASTA Nucleotide Changer

Converts FASTA sequences from/to RNA/DNA

## FASTQ Quality Filter

Filters sequences based on quality

## FASTQ Quality Trimmer

Trims (cuts) sequences based on quality

## FASTQ Masker

Masks nucleotides with 'N' (or other character) based on quality



# The FASTX toolkit

2 run modes:

Through command-line interface on your computer

- No queues or server downtime to worry about
- No large data transfers

Through graphical interface on the GALAXY web portal

- Avoid “complicated” command structure
- No need for plotting software (e.g. gnuplot)
- Potentially large processing power through the cloud.


# GALAXY

<https://usegalaxy.org/>

- Easy to make and share customized workflows for data processing
- Embeds other software within the site (e.g. fastx toolkit, SAMtools, Picard, GATK)
- Now available on the Amazon cloud for \$
- However, no software for *de novo* assembly / clustering available by default
- And you still need to understand what the parameters *mean*.

# GALAXY

<https://usegalaxy.org/u/pierredewit/w/test>





# Data processing steps for the next four days

